

MEDFIRST

CRICOTHYROTOMY KIT

OVERVIEW OF OUR EMERGENCY AIRWAY KIT:

Our CE-certified kit is specifically designed for use in emergency situations, particularly during CICO (Can't Intubate, Can't Oxygenate), where conventional airway access is restricted.



This set is intended for cases requiring an emergency airway at the front of the neck, enabling the creation of an opening in the trachea through cricothyrotomy, ensuring immediate patient ventilation and oxygenation.

This CE-certified cricothyrotomy kit strictly complies with the 2015 guidelines of the DAS (Difficult Airway Society)

ADVANTAGES OF THE KIT

Speed and Efficiency: The kit enables extremely rapid emergency intervention, helping to restore the patient's breathing in just minutes.

Maximum Safety: Each component is carefully selected to ensure a safe and effective cricothyrotomy procedure, minimizing the risk of complications.

Ease of Use: Clear graphic instructions guide you through each step of the procedure.

Sterile Reference Material: Each kit is delivered in a single-use rigid package to maintain sterility and protect the components from damage or deformation.

KEY COMPONENTS OF THE CRICOTHYROTOMY KIT

- **10mm Scalpel:** Provides surgical precision for a clean incision, allowing quick access to the cricothyroid membrane.
- **Introducer:** Functions as a bougie, guiding the tracheostomy tube into the trachea after incision.
- **Tracheostomy Tube:** 6.0mm diameter with balloon, ensuring effective ventilation through a secure fit, preventing air leaks.
- **10ml Syringe:** Used to inflate the balloon, ensuring precise tube positioning.
- **Cotton Tape:** Secures the tracheostomy tube in place, preventing unwanted movement and maintaining stable ventilation for the patient.



Find our product with additional details on its [dedicated page](#).

More informations, contact-us

MEDFIRST

8 esplanade Compans Caffarelli 31000 TOULOUSE, France - RCS Toulouse B 831156641



+33 5.32.63.00.14



ventes@medfirst.care



www.medfirst.care